

Protocol: pediatric_head_*Prematura vid fullgångstid Endast 32 Kanal spole

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	Head First	Imaging Mode	2D
Patient Position	Supine	Pulse Sequence	Spin Echo
Coil Configuration	Nova 32ch Brain for 750w	Imaging Options	Seq, EDR, Fast, SS, ART, ARC
Plane	3-PLANE		
Series Description	Localizer		
SCAN TIMING		SCAN RANGE	
TE	Minimum	FOV	24.0
Number of Echoes	1	Slice Thickness	10.0
TR	1000.0	Slice Spacing	6.0
Receiver Bandwidth	83.33		
IMAGE ENHANCE		ACQ TIMING	
Filter Choice	None	Freq	256
		Phase	128
		Freq DIR	Unswap
		# of Acq. Before Pause	0
		Phase FOV	0.90
		Auto Shim	Auto
		Phase Correction	No
GATING/TRIGGER		USER CVS	
Auto Trigger Type	Off	User CV2	240.00
		User CV13	1.00
FMRI		MULTI-PHASE	
PSD Trigger	Internal	Seperate Series	0
View Order	Bottom/Up	Mask Phase	0
# of Repetitions REST	0	Mask Pause	0
# of Repetitions ACTIVE	0	Acoustic Reduction Level	3
SAT		DIFFUSION	
Tag Type	None	Recon All Images	On
TRICKS		CONTRAST	
Pause On/Off	On	Contrast Yes/No	No
Auto Subtract	0		
Auto SCIC	Off		
OTHERS			
Protocol Notes	OBS!!! Ingen kilkudde under spolen.		

Protocol: pediatric_head_*Prematura vid fullgångstid Endast 32 Kanal spole

Sag T2 ART	PATIENT POSITION		IMAGING PARAMETERS	
	Patient Entry	Head First	Imaging Mode	2D
	Patient Position	Supine	Pulse Sequence	FSE-XL
	Coil Configuration	32 Ch Head Coil	Imaging Options	EDR, TRF, Fast, ART, ARC
	Plane	OBLIQUE	SCAN RANGE	
	Series Description	Sag T2 ART	FOV	20.0
	SCAN TIMING		Slice Thickness	3.0
	Flip Angle	140	Slice Spacing	0.3
	TE	105.0	ACQ TIMING	
	Number of Echoes	1	Freq	300
	TR	4283.0	Phase	240
	Echo Train Length	20	Freq DIR	S/I
	Receiver Bandwidth	31.25	Fat Shift DIR	Normal (S)
	IMAGE ENHANCE		NEX	1.00
	Filter Choice	A	# of Acq. Before Pause	0
	GATING/TRIGGER		Phase FOV	0.90
	Auto Trigger Type	Off	Auto Shim	Auto
	FMRI		Phase Correction	Yes
	PSD Trigger	Internal	USER CVS	
	View Order	Bottom/Up	User CV22	1.00
	# of Repetitions REST	0	TR Min	3500.0
	# of Repetitions ACTIVE	0	TR Max	6000.0
	SAT		MULTI-PHASE	
	Tag Type	None	Seperate Series	0
	TRICKS		Mask Phase	0
	Pause On/Off	On	Mask Pause	0
	Auto Subtract	0	Acoustic Reduction Level	1
	Auto SCIC	2	DIFFUSION	
			Recon All Images	On
			CONTRAST	
			Contrast Yes/No	No

Sag T2 ART

Protocol: pediatric_head_*Prematura vid fullgångstid Endast 32 Kanal spole

Ax T2 NEO ART	PATIENT POSITION		IMAGING PARAMETERS	
	Patient Entry	Head First	Imaging Mode	2D
	Patient Position	Supine	Pulse Sequence	FSE-XL
	Coil Configuration	32 Ch Head Coil	Imaging Options	EDR, TRF, Fast, ZIP512, FR, ART, ARC
	Plane	OBLIQUE	SCAN RANGE	
	Series Description	Ax T2 NEO ART	FOV	19.2
	SCAN TIMING		Slice Thickness	1.0
	Flip Angle	111	Slice Spacing	Intleave
	TE	147.0	ACQ TIMING	
	Number of Echoes	1	Freq	192
	TR	8332.0	Phase	192
	Echo Train Length	15	Freq DIR	A/P
	Receiver Bandwidth	25.00	Fat Shift DIR	Normal (A)
	IMAGE ENHANCE		NEX	1.00
	Filter Choice	None	# of Acq. Before Pause	0
	GATING/TRIGGER		Phase FOV	1.00
	Auto Trigger Type	Off	Auto Shim	Auto
	FMRI		Phase Correction	Yes
	PSD Trigger	Internal	USER CVS	
	View Order	Bottom/Up	User CV18	1.00
	# of Repetitions REST	0	User CV22	1.00
	# of Repetitions ACTIVE	0	TR Min	2500.0
	SAT		TR Max	11000.0
	Tag Type	None	MULTI-PHASE	
	TRICKS		Seperate Series	0
	Pause On/Off	On	Mask Phase	0
	Auto Subtract	0	Mask Pause	0
	Auto SCIC	2	Acoustic Reduction Level	1
			DIFFUSION	
			Recon All Images	On
			CONTRAST	
			Contrast Yes/No	No

Ax T2 NEO ART

Protocol: pediatric_head_*Prematura vid fullgångstid Endast 32 Kanal spole

Ax T1 FLAIR ART	PATIENT POSITION		IMAGING PARAMETERS	
	Patient Entry	Head First	Imaging Mode	2D
	Patient Position	Supine	Pulse Sequence	FSE-XL
	Coil Configuration	32 Ch Head Coil	Imaging Options	Seq, EDR, TRF, Fast, ART, T1flair
	Plane	OBLIQUE	SCAN RANGE	
	Series Description	Ax T1 FLAIR ART	FOV	20.0
	SCAN TIMING		Slice Thickness	3.0
	Flip Angle	111	Slice Spacing	0.3
	TE	24.0	ACQ TIMING	
	Number of Echoes	1	Freq	296
	TR	2992.5	Phase	234
	TI	872	Freq DIR	A/P
	Echo Train Length	6	Fat Shift DIR	Normal (A)
	Receiver Bandwidth	41.67	NEX	2.00
	IMAGE ENHANCE		Phase FOV	0.90
	Filter Choice	None	Auto Shim	Auto
	GATING/TRIGGER		Phase Correction	No
	Auto Trigger Type	Off	USER CVS	
	FMRI		User CV12	1.00
	PSD Trigger	Internal	User CV22	1.00
	View Order	Bottom/Up	MULTI-PHASE	
	# of Repetitions REST	0	Seperate Series	0
	# of Repetitions ACTIVE	0	Mask Phase	0
	SAT		Mask Pause	0
	Tag Type	None	Acoustic Reduction Level	1
	TRICKS		DIFFUSION	
	Pause On/Off	On	Recon All Images	On
	Auto Subtract	0	CONTRAST	
	Auto SCIC	2	Contrast Yes/No	No

Ax T1 FLAIR ART

Protocol: pediatric_head_*Prematura vid fullgångstid Endast 32 Kanal spole

Ax_3D_SWIEPI	PATIENT POSITION		IMAGING PARAMETERS	
	Patient Entry	Head First	Imaging Mode	3D
	Patient Position	Supine	Pulse Sequence	SPGR
	Coil Configuration	32 Ch Head Coil	Imaging Options	EDR, Fast, ARC
	Plane	OBLIQUE	PSD Name	ksep12
	Series Description	Ax_3D_SWIEPI	SCAN RANGE	
	SCAN TIMING		FOV	24.0
	Flip Angle	20	Slice Thickness	2.0
	TE	Min Full	ACQ TIMING	
	Number of Echoes	1	Freq	448
	TR	Minimum	Phase	448
	Number of Shots	32	Freq DIR	Unswap
	Receiver Bandwidth	83.33	NEX	1.00
	IMAGE ENHANCE		Phase FOV	1.00
	Filter Choice	None	Auto Shim	Auto
	GATING/TRIGGER		Phase Correction	No
	Auto Trigger Type	Off	USER CVS	
	FMRI		User CV1	1.75
	PSD Trigger	Internal	User CV2	1.00
	View Order	Bottom/Up	User CV7	6.00
	# of Repetitions REST	0	MULTI-PHASE	
	# of Repetitions ACTIVE	0	Seperate Series	0
	SAT		Mask Phase	0
	Tag Type	None	Mask Pause	0
	TRICKS		DIFFUSION	
	Pause On/Off	On	Recon All Images	On
	Auto Subtract	0	CONTRAST	
	Auto SCIC	2	Contrast Yes/No	No

Ax_3D_SWIEPI

Protocol: pediatric_head_*Prematura vid fullgångstid Endast 32 Kanal spole

Cor T2 PROP ART	PATIENT POSITION		IMAGING PARAMETERS	
	Patient Entry	Head First	Imaging Mode	2D
	Patient Position	Supine	Pulse Sequence	PROPELLER
	Coil Configuration	Nova 32ch Brain for 750w	Imaging Options	EDR, TRF, Fast, ART, ARC
	Plane	OBLIQUE	SCAN RANGE	
	Series Description	Cor T2 PROP ART	FOV	20.0
	SCAN TIMING		Slice Thickness	3.0
	Flip Angle	142	Slice Spacing	0.3
	Number of Echoes	1	ACQ TIMING	
	TR	5645.0	Freq	336
	Echo Train Length	30	Freq DIR	S/I
	Receiver Bandwidth	41.67	Fat Shift DIR	Normal (I)
	IMAGE ENHANCE		NEX	2.50
	Filter Choice	None	# of Acq. Before Pause	0
	GATING/TRIGGER		Auto Shim	Auto
	Auto Trigger Type	Off	Phase Correction	No
	FMRI		USER CVS	
	PSD Trigger	Internal	User CV5	1.00
	View Order	Bottom/Up	TR Min	3000.0
	# of Repetitions REST	0	TR Max	10000.0
	# of Repetitions ACTIVE	0	MULTI-PHASE	
	SAT		Seperate Series	0
	SAT Location	I	Mask Phase	0
	Tag Type	None	Mask Pause	0
	TRICKS		Acoustic Reduction Level	1
	Pause On/Off	On	DIFFUSION	
	Auto Subtract	0	Recon All Images	On
	Auto SCIC	2	CONTRAST	
			Contrast Yes/No	No

Cor T2 PROP ART

Protocol: pediatric_head_*Prematura vid fullgångstid Endast 32 Kanal spole

Ax DTI 15R B800 ART	PATIENT POSITION		IMAGING PARAMETERS	
	Patient Entry	Head First	Imaging Mode	2D
	Patient Position	Supine	Pulse Sequence	Spin Echo
	Coil Configuration	Nova 32ch Brain for 750w	Imaging Options	EPI, DIFF, Asset, ART
	Plane	OBLIQUE	SCAN RANGE	
	Series Description	Ax DTI 15R B800 ART	FOV	22.0
	SCAN TIMING		Slice Thickness	3.6
	TE	Minimum	Slice Spacing	0.0
	Number of Echoes	1	ACQ TIMING	
	TR	8033.0	Freq	128
	Number of Shots	1	Phase	128
	IMAGE ENHANCE		Freq DIR	R/L
	Filter Choice	None	Phase FOV	1.00
	GATING/TRIGGER		Auto Shim	Auto
	Auto Trigger Type	Off	Phase Correction	Yes
	FMRI		USER CVS	
	PSD Trigger	Internal	User CV5	1.00
	View Order	Bottom/Up	User CV17	1.00
	# of Repetitions REST	0	TR Min	2000.0
	# of Repetitions ACTIVE	0	TR Max	15000.0
Ax DTI 15R B800 ART	SAT		MULTI-PHASE	
	Tag Type	None	Seperate Series	0
	TRICKS		Mask Phase	0
	Pause On/Off	On	Mask Pause	0
	Auto Subtract	0	Acoustic Reduction Level	3
	Auto SCIC	2	DIFFUSION	
			Optimized TE	Yes
			Diffusion Directions	Tensor
			Number of Diffusion Directions	15
			Number of T2 Images	2
			Dual Spin Echo	On
			Diffusion Tenser Processing Output	No Selection
			Recon All Images	Off
			CONTRAST	
			Contrast Yes/No	No

Protocol: pediatric_head_*Prematura vid fullgångstid Endast 32 Kanal spole

Ax DTI 2mm 64R ART	PATIENT POSITION		IMAGING PARAMETERS	
	Patient Entry	Head First	Imaging Mode	2D
	Patient Position	Supine	Pulse Sequence	Spin Echo
	Coil Configuration	32 Ch Head Coil	Imaging Options	EDR, Cla, EPI, DIFF, Asset, ART
	Plane	OBLIQUE	SCAN RANGE	
	Series Description	Ax DTI 2mm 64R ART	FOV	20.0
	SCAN TIMING		Slice Thickness	2.0
	TE	Minimum	Slice Spacing	0.0
	Number of Echoes	1	ACQ TIMING	
	TR	7021.0	Freq	100
	Number of Shots	1	Phase	100
	IMAGE ENHANCE		Freq DIR	R/L
	Filter Choice	None	Phase FOV	1.00
	GATING/TRIGGER		Auto Shim	Auto
	Auto Trigger Type	Off	Phase Correction	Yes
	FMRI		USER CVS	
	PSD Trigger	Internal	User CV5	1.00
	View Order	Bottom/Up	User CV17	1.00
	# of Repetitions REST	0	TR Min	2000.0
	# of Repetitions ACTIVE	0	TR Max	17000.0
	SAT		MULTI-PHASE	
	Tag Type	None	Seperate Series	0
	Fat/Water Saturation	Fat	Mask Phase	0
	TRICKS		Mask Pause	0
	Pause On/Off	On	Acoustic Reduction Level	3
	Auto Subtract	0	DIFFUSION	
	Auto SCIC	2	Optimized TE	Yes
			Diffusion Directions	Tensor
			Number of Diffusion Directions	64
			Number of T2 Images	6
			Dual Spin Echo	Off
			Diffusion Tenser Processing Output	No Selection
			Recon All Images	Off
			CONTRAST	
			Contrast Yes/No	No

Ax DTI 2mm 64R ART

Protocol: pediatric_head_*Prematura vid fullgångstid Endast 32 Kanal spole

--- SWI Rörig Pat ---	PATIENT POSITION		IMAGING PARAMETERS	
	Patient Entry	Head First	Imaging Mode	2D
	Patient Position	Supine	Pulse Sequence	Spin Echo
	Coil Configuration	32 Ch Head Coil	Imaging Options	Seq, EDR, Fast, SS, ART, ARC
	Plane	3-PLANE	SCAN RANGE	
	Series Description	--- SWI Rörig Pat ---	FOV	24.0
	SCAN TIMING		Slice Thickness	10.0
	TE	Minimum	Slice Spacing	6.0
	Number of Echoes	1	ACQ TIMING	
	TR	1000.0	Freq	256
	IMAGE ENHANCE		Phase	128
	Filter Choice	None	Freq DIR	Unswap
	GATING/TRIGGER		# of Acq. Before Pause	0
	Auto Trigger Type	Off	Phase FOV	0.90
	FMRI		Auto Shim	Auto
	PSD Trigger	Internal	Phase Correction	No
	View Order	Bottom/Up	USER CVS	
	# of Repetitions REST	0	User CV2	240.00
	# of Repetitions ACTIVE	0	User CV13	1.00
	SAT		MULTI-PHASE	
	Tag Type	None	Seperate Series	0
TRICKS			Mask Phase	0
Pause On/Off		On	Mask Pause	0
Auto Subtract		0	Acoustic Reduction Level	3
			DIFFUSION	
			Recon All Images	On
			CONTRAST	
			Contrast Yes/No	No

--- SWI Rörig Pat ---

Protocol: pediatric_head_*Prematura vid fullgångstid Endast 32 Kanal spole

Ax_2D_SWIEPI	PATIENT POSITION		IMAGING PARAMETERS	
	Patient Entry	Head First	Imaging Mode	2D
	Patient Position	Supine	Pulse Sequence	Gradient Echo
	Coil Configuration	32 Ch Head Coil	Imaging Options	EDR, EPI, ARC
	Plane	AXIAL	PSD Name	ksep12
	Series Description	Ax_2D_SWIEPI	SCAN RANGE	
	SCAN TIMING		FOV	22.0
	Flip Angle	70	Slice Thickness	3.0
	TE	Min Full	Slice Spacing	0.0
	Number of Echoes	1	ACQ TIMING	
	TR	Minimum	Freq	288
	Number of Shots	14	Phase	288
	Receiver Bandwidth	83.33	Freq DIR	R/L
	IMAGE ENHANCE		NEX	4.00
	Filter Choice	None	Phase FOV	1.00
	GATING/TRIGGER		Auto Shim	Auto
	Auto Trigger Type	Off	Phase Correction	No
	FMRI		USER CVS	
	PSD Trigger	Internal	User CV1	1.00
	View Order	Bottom/Up	User CV2	1.00
	# of Repetitions REST	0	User CV3	1.00
	# of Repetitions ACTIVE	0	User CV7	6.00
	SAT		MULTI-PHASE	
	Tag Type	None	Seperate Series	0
	TRICKS		Mask Phase	0
	Pause On/Off	On	Mask Pause	0
	Auto Subtract	0	DIFFUSION	
	Auto SCIC	2	Recon All Images	On
			CONTRAST	
			Contrast Yes/No	No

Ax_2D_SWIEPI

Protocol: pediatric_head_*Prematura vid fullgångstid Endast 32 Kanal spole

--Eventuellt--	PATIENT POSITION		IMAGING PARAMETERS	
	Patient Entry	Head First	Imaging Mode	2D
	Patient Position	Supine	Pulse Sequence	Gradient Echo
	Coil Configuration	GEM Head	Imaging Options	Seq, Fast
	Plane	3-PLANE	SCAN RANGE	
	Series Description	--Eventuellt--	FOV	24.0
	SCAN TIMING		Slice Thickness	10.0
	Number of Echoes	1	Slice Spacing	5.0
	IMAGE ENHANCE		ACQ TIMING	
	Filter Choice	None	Freq	256
	GATING/TRIGGER		Phase	128
	Auto Trigger Type	Off	Freq DIR	Unswap
	MULTI-PHASE		NEX	1.00
	Seperate Series	0	Phase FOV	1.00
	Mask Phase	0	Auto Shim	Auto
	Mask Pause	0	Phase Correction	No
	DIFFUSION		FMRI	
	Recon All Images	On	PSD Trigger	Internal
	CONTRAST		Slice Order	Interleaved
	Contrast Yes/No	No	View Order	Bottom/Up
			# of Repetitions REST	0
			# of Repetitions ACTIVE	0
			SAT	
			Tag Type	None
			TRICKS	
			Pause On/Off	On
			Auto Subtract	0
			Auto SCIC	Off

--Eventuellt--

Protocol: pediatric_head_*Prematura vid fullgångstid Endast 32 Kanal spole

Ax T1 FLAIR PROP ART	PATIENT POSITION		IMAGING PARAMETERS	
	Patient Entry	Head First	Imaging Mode	2D
	Patient Position	Supine	Pulse Sequence	PROPELLER
	Coil Configuration	Nova 32ch Brain for 750w	Imaging Options	EDR, TRF, Fast, ART, ARC, T1flair
	Plane	OBLIQUE	SCAN RANGE	
	Series Description	Ax T1 FLAIR PROP ART	FOV	22.0
	SCAN TIMING		Slice Thickness	3.0
	Flip Angle	111	Slice Spacing	0.3
	Number of Echoes	1	ACQ TIMING	
	TR	2802.0	Freq	280
	TI	1131	Freq DIR	A/P
	Echo Train Length	14	Fat Shift DIR	Normal (A)
	Receiver Bandwidth	41.67	NEX	3.50
	IMAGE ENHANCE		# of Acq. Before Pause	0
	Filter Choice	None	Auto Shim	Auto
	GATING/TRIGGER		Phase Correction	No
	Auto Trigger Type	Off	USER CVS	
	FMRI		User CV5	1.00
	PSD Trigger	Internal	User CV12	1.00
	View Order	Bottom/Up	MULTI-PHASE	
	# of Repetitions REST	0	Seperate Series	0
	# of Repetitions ACTIVE	0	Mask Phase	0
	SAT		Mask Pause	0
	Tag Type	None	Acoustic Reduction Level	1
	TRICKS		DIFFUSION	
	Pause On/Off	On	Recon All Images	On
	Auto Subtract	0	CONTRAST	
	Auto SCIC	2	Contrast Yes/No	No

Ax T1 FLAIR PROP ART

Protocol: pediatric_head_*Prematura vid fullgångstid Endast 32 Kanal spole

Ax T2 PROP ART	PATIENT POSITION		IMAGING PARAMETERS	
	Patient Entry	Head First	Imaging Mode	2D
	Patient Position	Supine	Pulse Sequence	PROPELLER
	Coil Configuration	32 Ch Head Coil	Imaging Options	EDR, TRF, Fast, FR, ARC
	Plane	OBLIQUE	SCAN RANGE	
	Series Description	Ax T2 PROP ART	FOV	20.0
	SCAN TIMING		Slice Thickness	3.0
	Flip Angle	111	Slice Spacing	0.3
	Number of Echoes	1	ACQ TIMING	
	TR	8852.0	Freq	360
	Echo Train Length	31	Freq DIR	A/P
	Receiver Bandwidth	41.67	Fat Shift DIR	Normal (A)
	IMAGE ENHANCE		NEX	2.30
	Filter Choice	None	Auto Shim	On
	GATING/TRIGGER		Phase Correction	No
	Auto Trigger Type	Off	USER CVS	
	FMRI		User CV5	1.00
	PSD Trigger	Internal	TR Min	400.0
	View Order	Bottom/Up	TR Max	11000.0
	# of Repetitions REST	0	MULTI-PHASE	
	# of Repetitions ACTIVE	0	Seperate Series	0
	SAT		Mask Phase	0
	Tag Type	None	Mask Pause	0
	TRICKS		DIFFUSION	
	Pause On/Off	On	Recon All Images	On
	Auto Subtract	0	CONTRAST	
	Auto SCIC	2	Contrast Yes/No	No

Ax T2 PROP ART